


MS Spectrum



Sample Information

Dissolution method	:15% ACN+85% H ₂ O	Interface	:ESI	Prerod Bias	:+4.5kv
Date Acquired	:2023/12/22 13:20:49	Nebulizing Gas Flow	:1.5L0/min	Detector	:-0.2kv
Injection Volume	:0.2ul	CDL Temp	:250°C	T.Flow	:0.2ml/min
Block Temp	:200	CDL Volt	:0v	B.conc	:50% H ₂ O/50% MEOH
Name	:Semaglutide				
Sequence	:His-Aib-Glu-Gly-Thr-Phe-Thr-Ser-Asp-Val-Ser-Ser-Tyr-Leu-Glu-Gly-Gln-Ala-Ala-Lys(AEEA-AEEA-γ-Glu-17-carboxyheptadecanoyl)-Glu-Phe-Ile-Ala-Trp-Leu-Val-Arg-Gly-Arg-Gly-OH				
Theoretical	:4113.64				
Observed	:4113.75				

邮箱：slaes@motifbiotech.com

网址：<https://www.motifbiotech.com>

地址：苏州市高新区浒墅关镇浒青路 111 号 B04