


MS Spectrum

Positive
Intensity



Sample Information

Dissolution method	:15% ACN+85% H ₂ O	Interface	:ESI	Prerod Bias	:+4.5kv
Date Acquired	:2024/10/21 14:15:01	Nebulizing Gas Flow	:1.50L/min	Detector	:-0.2kv
Injection Volume	:0.2ul	CDL Temp	:250°C	T.Flow	:0.2ml/min
Block Temp	:200	CDL Volt	:0v	B.conc	:50% H ₂ O/50% MEOH

Name	:TfR Targeting Peptide
Sequence	:THRPPMWSPVWP
Theoretical	:1490.73
Observed	:1490.70