

CERTIFICATE OF ANALYSIS

| | |
|------------------------------|----------------------------|
| Product Name | Acetyl tetrapeptide-3 |
| Sequence | Ac-KGHK-NH2 (Acetate Salt) |
| Dissolution condition | 100%H2O |
| Length | 4AA |
| Modification | acetate salt |
| Molecular Weight (MW) | 509.61 |
| Storage | -20°C |

| Test Items | Specifications | Results |
|----------------------------|---------------------------------------|----------------|
| MW by MS | 509.70 | Conforms |
| Purity by HPLC | >95% | 97.804% |
| Peptide content | N/A | N/A |
| Moisture content | N/A | N/A |
| Acetic acid content | N/A | N/A |
| Appearance | White to off-white lyophilized powder | Conforms |
| Quantity | 2g | 50.0mg*40Tube |

Certified by:

Quality Assurance Department



Note: this product is intended for research use only; not for diagnostic or human use.

Mail : marketing@motifbiotech.com

Web : <https://www.motifbiotech.com>

Add : B04, Building 17, No. 122, Yongan Rd, Suzhou National Hi-Tech District, PRC