

CERTIFICATE OF ANALYSIS

Product Name	Acetyl Tetrapeptide-9
Sequence	Ac-QDVH
Dissolution condition	100%H ₂ O
Length	4AA
Modification	acetate salt
Molecular Weight (MW)	539.55
Storage	-20°C

Test Items	Specifications	Results
MW by MS	540.35(Peptide)	Conforms
Purity by HPLC	>98%	98.99%
Peptide content	N/A	83.39%
Moisture content	7.0%	2.24%
Acetic acid content	15.0%	8.55%
Appearance	White to off-white powd	Conforms
Quantity	10g	10.0g



Certified by:
Quality Assurance Department

Note: this product is intended for research use only; not for diagnostic or human use.

Mail : marketing@motifbiotech.com

Web : <https://www.motifbiotech.com>

Add : B04, Building 17, No. 122, Yongan Rd, Suzhou National Hi-Tech District, PRC