


MS Spectrum



Sample Information

Dissolution method :15% ACN+85% H<sub>2</sub>O  
Date Acquired :2024/04/23 16:40:34  
Injection Volume :0.2ul  
Block Temp :200

Interface :ESI  
Nebulizing Gas Flow :1.50L/min  
CDL Temp :250°C  
CDL Volt :0v

Prerod Bias :+4.5kv  
Detector :-0.2kv  
T.Flow :0.2ml/min  
B.conc :50% H<sub>2</sub>O/50% MEOH

Name :Biotin Tripeptide-1  
Sequence :Biotin-GHK  
Theoretical :566.69  
Observed :566.40

Email : marketing@motifbiotech.com

Web : <https://www.motifbiotech.com>

Add : B04, Building 17, No. 122, Yongan Rd, Suzhou National Hi-Tech District, PRC