

CERTIFICATE OF ANALYSIS

Product Name	Decapeptide-12
Sequence	YRSRKYSSWY (Acetate Salt)
Dissolution condition	100%H ₂ O
Length	10AA
Modification	acetate salt
Molecular Weight (MW)	1395.52
Storage	-20°C

Test Items	Specifications	Results
MW by MS	1395.30	Conforms
Purity by HPLC	>95%	98.926%
Peptide content	N/A	N/A
Moisture content	N/A	N/A
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	1g	50.0mg*20Tube

Certified by:

Quality Assurance Department

**Note: this product is intended for research use only; not for diagnostic or human use.**Mail : marketing@motifbiotech.comWeb : <https://www.motifbiotech.com>

Add : B04, Building 17, No. 122, Yongan Rd, Suzhou National Hi-Tech District, PRC