

CERTIFICATE OF ANALYSIS

Product Name	Dipeptide-6
Sequence	Pro-Hyp
Dissolution condition	100%H ₂ O
Length	2AA
Modification	Hyp
Molecular Weight (MW)	228.26
Storage	-20°C

Test Items	Specifications	Results
MW by MS	228.10	Conforms
Purity by HPLC	>98%	98.355%
Peptide Content	N/A	N/A
Moisture content	N/A	N/A
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	100mg	100.0mg

Certified by:

Quality Assurance Department

**Note: this product is intended for research use only; not for diagnostic or human use.**Mail : marketing@motifbiotech.comWeb : <https://www.motifbiotech.com>

Add : B04, Building 17, No. 122, Yongan Rd, Suzhou National Hi-Tech District, PRC