

CERTIFICATE OF ANALYSIS

Product Name	Hexapeptide-2
Sequence	H-{D-Trp}-AW-{D-Phe}-K-NH ₂
Dissolution condition	100%H ₂ O
Length	6AA
Modification	sodium salt
Molecular Weight (MW)	873.01
Storage	-20°C

Test Items	Specifications	Results
MW by MS	873.92(Peptide)	Conforms
Purity by HPLC	>98%	99.70%
Peptide content	N/A	N/A
Moisture content	7.0%	1.29%
Acetic acid content	N/A	N/A
Appearance	White to off-white powd	Conforms
Quantity	2g	2.0g



Certified by:
Quality Assurance Department

Note: this product is intended for research use only; not for diagnostic or human use.

Mail : marketing@motifbiotech.com

Web : <https://www.motifbiotech.com>

Add : B04, Building 17, No. 122, Yong'an Rd, Suzhou National Hi-Tech District, PRC