

CERTIFICATE OF ANALYSIS

Product Name	Palmitoyl Tripeptide-1
Sequence	Pal-GHK (HCl Salt)
Dissolution condition	50%ACN+50%H2O
Length	3AA
Modification	HCl Salt
Molecular Weight (MW)	578.79
Storage	-20°C

Test Items	Specifications	Results
MW by MS	578.80	Conforms
Purity by HPLC	>98%	98.141%
Peptide content	N/A	N/A
Moisture content	N/A	N/A
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	1g	1.0g

Certified by:

Quality Assurance Department

**Note: this product is intended for research use only; not for diagnostic or human use.**Mail : marketing@motifbiotech.comWeb : <https://www.motifbiotech.com>

Add : B04, Building 17, No. 122, Yong'an Rd, Suzhou National Hi-Tech District, PRC