


MS Spectrum



Sample Information

Dissolution method	:5% HAC+8% ACN+87% H2O	Interface	:ESI	Prerod Bias	:+4.5kv
Date Acquired	:2023/11/18 16:25:57	Nebulizing Gas Flow	:1.50L/min	Detector	:-0.2kv
Injection Volume	:0.2ul	CDL Temp	:250C	T.Flow	:0.2ml/min
Block Temp	:200	CDL Volt	:0v	B.conc	:50% H2O/50% MEOH

Name	:Palmitoyl Tripeptide-1
Sequence	:Pal-GHK (HCl Salt)
Theoretical	:578.79
Observed	:578.80