

CERTIFICATE OF ANALYSIS

Product Name	Pentapeptide-3
Sequence	GPRPA-NH2 (Acetate Salt)
Dissolution condition	100%H2O
Length	5AA
Modification	Acetate salt
Molecular Weight (MW)	495.56
Storage	-20°C

Test Items	Specifications	Results
MW by MS	495.60	Conforms
Purity by HPLC	>95%	97.649%
Peptide Content	N/A	N/A
Moisture content	N/A	N/A
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	40mg	5.0mg*8Tube

Certified by:

Quality Assurance Department

**Note: this product is intended for research use only; not for diagnostic or human use.**Email : marketing@motifbiotech.comWeb : <http://www.motifbiotech.com>

Add : B04, Building 17, No. 122, Yongan Rd, Suzhou National Hi-Tech District, PRC