

CERTIFICATE OF ANALYSIS

Product Name	Tripeptide-41
Sequence	FWY (Acetate Salt)
Dissolution condition	5%HCOOH+15%ACN+80%H ₂ O
Length	3AA
Modification	acetate salt
Molecular Weight (MW)	514.58
Storage	-20°C

Test Items	Specifications	Results
MW by MS	514.40	Conforms
Purity by HPLC	>95%	99.796%
Peptide content	N/A	N/A
Moisture content	N/A	N/A
Acetic acid content	N/A	N/A
Appearance	White to off-white lyophilized powder	Conforms
Quantity	1g	20.0mg*50Tube

Certified by:

Quality Assurance Department

**Note: this product is intended for research use only; not for diagnostic or human use.**

Mail : marketing@motifbiotech.com

Web : <https://www.motifbiotech.com>

Add : B04, Building 17, No. 122, Yong'an Rd, Suzhou National Hi-Tech District, PRC