


MS Spectrum



Sample Information

| | | | | | |
|--------------------|-------------------------------|---------------------|------------|-------------|--------------------------------|
| Dissolution method | :15% ACN+85% H ₂ O | Interface | :ESI | Prerod Bias | :+4.5kv |
| Date Acquired | :2024/09/26 10:39:30 | Nebulizing Gas Flow | :1.50L/min | Detector | :-0.2kv |
| Injection Volume | :0.2ul | CDL Temp | :250°C | T.Flow | :0.2ml/min |
| Block Temp | :200 | CDL Volt | :0v | B.conc | :50% H ₂ O/50% MEOH |

| | |
|-------------|-------------------------|
| Name | :Palmitoyl tripeptide-5 |
| Sequence | :Pal-KVK |
| Theoretical | :611.90 |
| Observed | :611.80 |