




Sample Information

Dissolution method	:15%ACN+85%H ₂ O	Interface	:ESI	Prerod Bias	:+4.5kv
Date Acquired	:2025/03/12 14:54:44	Nebulizing Gas Flow	:1.50L/min	Detector	:-0.2kv
Injection Volume	:0.2ul	CDL Temp	:250°C	T.Flow	:0.2ml/min
Block Temp	:200	CDL Volt	:0v	B.conc	:50%H ₂ O/50%MEOH
Name	:卡格列肽				
Sequence	:{Eicosanedioic acid-γ-Glu}- KC*NTATC*ATQRLAEFL*HSSNFGPILPPTNVGSNTP-NH ₂ (Acetate Salt)				
Theoretical	:4409.01				
Observed	:4409.40				

邮箱: marketing@motifbiotech.com

网址: <https://www.motifbiotech.com>

地址: 苏州市高新区浒墅关镇浒青路 111 号 B04